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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

## Brace for acromioclavicular dislocation

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

### SELECTION/SIZE

Code	REF.8501
Size	universal
Colour	beige

### MAINTENANCE

- ⌘ Do not bleach
- ⌘ No chemical cleaning
- ⌘ Do not iron
- ⌘ Do not tumble-dry
- ⌘ Washing instructions: hand wash at max. 30°C with neutral soap.

Do not dispose of the product or any of its components into the environment.

### INDICATIONS

- Stabilises the shoulder and arm in cases of dislocation and subluxation of the acromioclavicular joint (grade I and II)

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

- Soft felt support belts
- Straps with locking hooks and buckles
- Cotton blend fabric structure

### PUTTING ON THE APPLIANCE

- 1 The package contains the arm bag (C) and the two padded straps: the shoulder strap (shorter) (A) and the shoulder strap (longer) (B).
- 2 Position the shoulder strap (A) over the affected shoulder, as shown in the figure.
- 3 Hook the shoulder strap back to the shoulder strap (D). Pass the shoulder strap under the armpit, making sure that the hook (E) is facing the low and hook at the front (F).
- 4 Gently house the flexed arm inside the canvas bag (C).
- 5 Attach the bag buckles to the hook (E) on the shoulder strap and to the hooks (G) and (H) on the shoulder strap.
- 6 Tighten the tie rods as needed. Raise your elbow as you increase the tension on the front and back buckles.

**Note:** The under-shoulder tie-rod is equipped with two hooks (E-I): when used on the injured right arm (see figure) only the hook (E) is used, when used on the left injured arm only the hook (I).

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Velcro® est une marque enregistrée de Velcro Industries B.V.