

Nome Prodotto

Ruote da Transito.

| Codice | Ø (cm) |
|---------------|----------|
| 01.04.12QRPUS | 30 (12") |
| 01.04.12QRPUA | 30 (12") |

**Destinazione d'uso**

Le ruote da transito vengono utilizzate sulle carrozzine con ruote posteriori ad estrazione rapida (Quick Release).
La loro funzione è di sostituire le classiche ruote da 60 cm (24") con ruote ad ingombro ridotto (Ø 30 cm, 12").

| Compatibilità | |
|--|--|
| 01.04.12QRPUS | Carrozzine Viky Exchange • Tutte le Misure |
| 01.04.12QRPUA | Carrozzine Viky Comfy Light • Tutte le Misure |
| Materiali | |
| 01.04.12QRPUS | PU - Acciaio |
| 01.04.12QRPUA | PU - Alluminio |
| Manutenzione e Pulizia | |
| Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica | |
| Conservazione e Smaltimento | |
| Conservare al riparo da fiamme libere e fonti di calore | |
| Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata | |
| Confezione e Nr. Articoli | |
| Scatola | 2 Pezzi |
| Garanzia | |
| Europea | 24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili |
| Riferimenti del Prodotto | |
| Fabbricante | Demarta Virginio S.a.S. |
| Classe di Dispositivo 93/42 | Classe I |
| Codice CND | Y1224 |

**Product Name**

Transit Wheels.



| Code | Ø (cm) |
|---------------|----------|
| 01.04.12QRPUS | 30 (12") |
| 01.04.12QRPUA | 30 (12") |

Use

The transit wheels are used on wheelchairs with quick-release rear wheels.

Their function is to replace the normal 60 cm (24 ") wheels with smaller one (Ø 30 cm, 12 ").

| Compatibility | |
|--|--|
| 01.04.12QRPUS | Viky Exchange Wheelchair • All Dimensions |
| 01.04.12QRPUA | Viky Comfy Light Wheelchair • All Dimensions |
| Materials | |
| 01.04.ATS | PU - Steel |
| 01.04.ATA | PU - Aluminum |
| Maintenance and Cleaning | |
| Disinfection is possible with 5% diluted disinfectants even with alcohol base | |
| Storage and Disposal | |
| Keep away from heat, sparks, open flame or other sources of ignition | |
| Dispose of the device in compliance with current regulations on environmental protection and separate waste collection | |
| Package and Nr.Items | |
| Box | 2 Pieces |
| Warranty | |
| European | 24 months for manufacturing defects with the exception of consumable parts |
| Product References | |
| Manufacturer | Demarta Virginio S.a.S. |
| Directive 93/42 | Classe I |
| CND Code | Y1224 |