

**Nome Prodotto**

Schienale Tensionabile.

Codice	Larghezza (cm)
01.04.VK513T	41
01.04.VK514T	43
01.04.VK515T	46
01.04.VK516T	48

**Destinazione d'uso**

Lo Schienale Tensionabile segue la forma della colonna dell'utilizzatore, è personalizzabile in base alle necessità. E' molto utile per accogliere e sostenere deformità tipo cifosi o inclinazioni pelviche.

Compatibilità	
Carrozine Viky Comfy Light	
Materiali	
Schienale Tensionabile	Poliestere
Manutenzione e Pulizia	
Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica	
Conservazione e Smaltimento	
Si consiglia l'uso di un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.a.S.
Classe di Dispositivo 93/42	Classe I
Codice ISO	12.22.91.709
Codice CND	Y1224



**Product Name**

Tensionable Backrest.

Code	Width (cm)
01.04.VK513T	41
01.04.VK514T	43
01.04.VK515T	46
01.04.VK516T	48



**Use**

The Tensionable Backrest follows the shape of the user's spine, and can be customized according to need. It is very useful to accommodate and support kyphosis deformities or pelvic inclinations.

Compatibility	
Viky Comfy Light model	
Materials	
Tensionable Back	Polyester
Maintenance and Cleaning	
Disinfection is possible with 5% diluted disinfectants even with alcohol base	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.a.S.
Directive 93/42	Classe I
ISO Code	12.22.91.709
CND Code	Y1224