

**Sedia top class
disegnata dal vento**



EUREKA SC

kg
12




EUREKA SUPER

kg
12,5



Schienale estraibile (optional)
Removable backrest (optional)



SURACE INNOVATION Sedie leggere Light weight wheelchairs

EUREKA
Sedia top class leggerissima

Standard: Telaio modulare integralmente in alluminio con doppia crociera. Innovativo sistema facilitato di apertura e chiusura. Possibilità di regolazione del sedile in altezza. Braccioli desk ribaltabili. Astine porta pedane girevoli, estraibili, senza ingombro sul montante anteriore. Ruote in alluminio ad estrazione rapida. Ruote ascensore ed antiribaltamento integrate.

Standard sul modello EUREKA SC: Fianco giro ruota in ABS rinforzato. Ruote anteriori con forcelle ad angolo variabile.

Standard sul modello EUREKA SUPER: Ruote sport 24x1 e cannotti forcelle ad angolo variabile.

Optional: Braccioli regolabili in altezza ed estraibili. Cannotti forcelle anteriori ad angolo variabile. Appoggiambraccioli regolabili in altezza, girevoli ed estraibili. Schienale tensionabile. Schienale estraibile. Schienale regolabile in inclinazione. Monoguida. Pedana unica. Tavolino sagomato. Cintura di sicurezza. Gomme sport 24x1.

EUREKA
The lightest top class chair

Standard: All aluminium, modular frame with double cross. Innovative, easy opening and closure system. Possibility of height adjustment of seat. Desk and tipping armrests. Swinging and detachable footrests without any obstacle on front support bar. Quickly extractable and aluminium big wheels. Integrated wheels for lift and anti-overturn.

Standard on model EUREKA SC: Side-armrest contouring wheels in special ABS. Front wheels with angle adjustable forks.

Standard on model EUREKA SUPER: 24x1 sport wheels and angle adjustable forks.

Optional: Adjustable in height and removable armrests. Angle adjustable, front forks. Adjustable in height, swinging and detachable legrests. Tension adjustable backrest. Removable backrest. Reclining backrest. One hand driving system. Unique footrest. Contoured table. Security belt. 24x1 sport tyres.

ISO Class: 12.21.06.060 o: 12.21.06.039 + 12.24.03.121 + 12.24.21.106 + 12.24.03.118 + 12.24.21.109 + 12.24.21.112 + 12.24.06.172

Codice LEA 12.22.03.012

