

## Per la vita dinamica



kg 6,5



+

kg 1,5



+

kg 1,5



+

kg 3,5



=

+ kg 1  
Ruote antiribaltamento



**S13**

### Sedia super compact leggera

**Standard:** Telaio in alluminio ed acciaio speciale. Doppia crociera. Braccioli desk regolabili in altezza. Fianchi estraibili. Pedane girevoli ed estraibili. Ruote ad estrazione rapida. Rotelline per ascensore ed antiribaltamento. Regolazione sedile in altezza. Catarifrangenti.

**Optional:** Appoggiamambe regolabili girevoli ed estraibili. Schienale tensionabile.



**S13**

### Light and super compact chair

**Standard:** Frame made of aluminium and special steel. Double cross. Adjustable and removable desk arms. Swinging and detachable footrests. Quickly extractable big wheels. Small wheels for lift and anti-overturn. Seat rest adjustable in height. Reflectors.

**Optional:** Adjustable in height and removable legrests. Tension adjustable backrest.

Codice LEA 12.22.03.012

ISO Class 12.21.06.039+12.24.21.106+


12.24.03.121+12.24.03.118




cm Ø 60



cm Ø 18



L 42



L 61



L 75



L 96



H 50



Standard



Portata

130 kg