

Leggera e indistruttibile Light and unbreakable



Braccioli regolabili.
Kit speed (optional)
Adjustable armrests.
Kit speed (optional)



■ SUPERITALA Sedia leggera

Standard: Telaio in alluminio ed acciaio speciale ad alta resistenza. Ruote pneumatiche ad estrazione rapida con ruote ascensore. Monocrociera rinforzata. Braccioli desk girevoli. Pedane regolabili, girevoli ed estraibili. Regolazione del sedile in altezza.

Optional: Kit speed (passo accorciato e rotelline antiribaltamento). Braccioli regolabili in altezza. Catarifrangenti. Cintura di sicurezza.

🇬🇧 SUPERITALA Light chair

Standard: Aluminium and special highstrength steel frame. Pneumatic extractable big wheels and wheels for lift. Single reinforced cross. Swinging and desk armrests. Swinging, adjustable and removable footrests. Adjustable in height seat.

Optional: Speed kit (shortened step with anti-overturn small wheels). Adjustable in height armrests. Reflectors. Security belt.

Codice LEA 12.22.03.006 (ricond.)

ISO Class 12.21.06.039 + 12.24.03.121 + 12.24.21.106




Telaio
stabilizzante standard




cm Ø 60



cm Ø 18



L 43



L 61



h 51



L 75 - 100



L 69 - 94
con Kit Speed



Portata
130 kg