

Sedia in alluminio Aluminium chair



Braccioli regolabili in altezza (optional)
Adjustable in height armrests (optional)



Appoggiamambe regolabili (optional)
Adjustable in height legrests (optional)

EUREKINA
Sedia da transito leggerissima

Standard: Telaio modulare integralmente in alluminio. Braccioli desk ribaltabili. Doppia crociera.

Optional: Appoggiamambe regolabili, girevoli ed estraibili. Braccioli regolabili in altezza ed estraibili. Tavolino sagomato. Tasca allo schienale. Cintura di sicurezza.

EUREKINA
The lightest transit chair

Standard: All aluminium, modular frame. Desk and tipping armrests. Double cross.

Optional: Adjustable in height, swinging and detachable legrests. Adjustable in height and removable armrests. Contoured table. Pocket at the back. Security belt.

kg 10,7 = Kg 9,6 + Kg 1,1

Codice LEA 12.22.18.003

ISO Class: 12.21.06.045 + 12.24.03.121 + 12.24.03.118 + 12.24.06.172

cm Ø 30	cm Ø 18	L 44 - 54	L 68 - 95	h 51	h 36 - 44	Standard	Optional	Optional	Optional	Portata 130 kg